



Table of Contents

| | | |
|-------|---|----|
| 1 | Principle and scope..... | 3 |
| 2 | Validation..... | 3 |
| 3 | Reagents, chemicals and reference materials..... | 4 |
| 3.1 | Reagents and Chemicals..... | 4 |
| 3.2 | Reference Materials..... | 4 |
| 4 | Equipment..... | 5 |
| 5 | Solutions..... | 6 |
| 5.1 | Dilution solvent..... | 6 |
| 5.2 | Preparation of reference stock solutions..... | 6 |
| 5.3 | Preparation of comparison solutions and LOD solution..... | 6 |
| 5.4 | Preparation of blank solution..... | 6 |
| 5.5 | Preparation of test solutions..... | 7 |
| 6 | Experimental..... | 7 |
| 6.1 | HPLC..... | 7 |
| 6.2 | Single quad mass spectrometer..... | 8 |
| 6.3 | Order of injections..... | 8 |
| 7 | System Suitability Test..... | 9 |
| 7.1 | System Suitability Tests..... | 9 |
| 7.1.1 | SST - LOQ - Signal to noise ratio..... | 9 |
| 7.1.2 | SST - Performance control..... | 9 |
| 8 | Evaluation and reporting..... | 9 |
| 9 | Chromatograms..... | 10 |
| 9.1 | Blank solution..... | 10 |
| 9.2 | Comparison solution 1 (CS1) / LOQ / 0.5 ppm..... | 10 |
| 9.3 | Test solution (epoxy resin)..... | 10 |
| 10 | Abbreviations..... | 11 |
| 11 | Changes to previous edition..... | 11 |

1 Principle and scope

This method is used to determine Bisphenol A in different BPA-epoxy resins.

The determination of Bisphenol A is carried out by liquid chromatography-mass spectrometry (LC/MS). The specific mass of Bisphenol A is measured in the selected ion monitoring mode (SIM-mode).



MW 228,28

The method is developed as a quantitative test method:

| Test | Range [ppm, mg/kg] ^[A] | Limit of quantification (LOQ) [ppm, mg/kg] ^[A] | Limit of detection (LOD) [ppm, mg/kg] ^[A] |
|-------------|--------------------------------------|--|---|
| Bisphenol A | 0.5 - 10 | 0.5 | 0.25 |

^[A] Related to 75 mg test sample and a final volume of 10 mL

2 Validation

This quantitative test method was validated according to the ICH Guideline Q2 (R1). It includes:

- Accuracy
- Linearity
- Specificity
- Precision (repeatability)
- Intermediate precision
- Limit of quantitation (LOQ)
- Limit of detection (LOD)
- Stability of solutions

**BPA-epoxy resins:
Determination of Bisphenol A by LC-MS**

C.52.S1593_02

3 Reagents, chemicals and reference materials
3.1 Reagents and Chemicals

The use of chemicals and reagents of equivalent or higher quality, also delivered by differing than mentioned suppliers, is allowed.

| Reagent/Chemical | Quality | Supplier | Product Number |
|-----------------------------|--|--------------------|----------------|
| Ammonium acetate | for HPLC, $\geq 99.0\%$ | Fluka | 17836 |
| Acetonitrile | CHROMASOLV®, gradient grade, for HPLC, $\geq 99.9\%$ | Sigma Aldrich | 34851 |
| Ammonium hydroxide solution | $\geq 25\%$ in H ₂ O, eluent additive for LC-MS | Fluka | 44273 |
| Water deionized | 18.2 MΩcm | Millipore, Milli Q | n/a |

3.2 Reference Materials

| Reagent/Chemical | Quality | Supplier | Part. no. |
|------------------|-------------|---------------|-----------|
| Bisphenol A | $\geq 99\%$ | Sigma-Aldrich | 239658 |

| CAS No. | Toxicity | Synonyms |
|---------|--|---|
| 80-05-7 | Solvias category 2 (H317, H318, H335, H361, H411) | 2,2-Bis(4-hydroxyphenyl)propane, 4,4'-Isopropylidenediphenol |

4 Equipment

Use the following equipment or equipment of equivalent performance.

| Instrument / Accessories | Manufacturer / Provider | Model |
|---------------------------------------|-------------------------|---|
| High performance liquid chromatograph | Agilent | 1200 Series |
| Mass spectrometer (MS) | Agilent | MSD6140 |
| MS Electro Spray Ion Source | Agilent | G1948A |
| Chromatography data system (CDS) | Agilent | LC/MSD ChemStation, Rev. B.04.03-SP2 |
| Eppendorf pipette | Eppendorf | 10 µL – 100 µL |
| Eppendorf pipette | Eppendorf | 100 µL – 1000 µL |
| Eppendorf pipette | Eppendorf | 500 µL – 5000 µL |
| pH meter | Metrohm | pH meter 780 |
| Analytical balance | Mettler Toledo | XP205 |
| Precision balance | Mettler Toledo | XP2002S |
| Standard crimp vial 2 mL, amber glass | Infochroma | 8C02-CV(A) |
| Alu crimp cap with rubber/PTFE septa | Infochroma | G003-AC-RU/TE |
| UptiDisc Syringe filter | Interchim | PVDF or PTFE / 25 mm / 0.45 µm |

5 Solutions

Weights and volumes may be adapted as long as the concentrations remain unchanged.

5.1 Dilution solvent

Pipette 50 mL of acetonitrile in a 100 mL volumetric flask and dilute to volume with water.

5.2 Preparation of reference stock solutions

Reference stock solution 1 (RSS1)

Weigh 48 mg to 52 mg of Bisphenol A (accuracy: 0.01 mg) in a 50 mL volumetric flask. Dissolve and dilute to volume with acetonitrile.

(Concentration of Bisphenol A: 1000 µg/mL)

Reference stock solution 2 (RSS2)

Pipette 1000 µL of reference stock solution 1 (RSS1) into a 10 mL volumetric flask and dilute to volume with acetonitrile.

(Concentration of Bisphenol A: 100 µg/mL)

Reference stock solution 3 (RSS3)

Pipette 150 µL of reference stock solution 2 (RSS2) into a 10 mL volumetric flask and dilute to volume with acetonitrile.

(Concentration of Bisphenol A: 1.5 µg/mL)

5.3 Preparation of comparison solutions and LOD solution

Pipette the amounts of reference stock solution 3 (RSS3) in a 10 ml volumetric flask according to the table below and dilute to volume with dilution solvent.

| Comparison solution | Addition of RSS3 [µL] | Conc. of Bisphenol A in solution [µg/mL] | Bisphenol A [ppm, mg/kg] ^[A] |
|---------------------|-----------------------|--|---|
| CS1 / LOQ | 25 | 0.0038 | 0.5 |
| CS2 | 50 | 0.0075 | 1.0 |
| CS3 | 125 | 0.0188 | 2.5 |
| CS4 | 250 | 0.0375 | 5.0 |
| CS5 | 375 | 0.0563 | 7.5 |
| CS6 | 500 | 0.0750 | 10.0 |

^[A] Related to 75 mg of test sample and a final volume of 10 mL

Transfer an aliquot of the clear solution in an amber vial and analyse by LC/MS.

The comparison solutions are stable for 3 days if stored at 20 °C and daylight in the autosampler tray.

5.4 Preparation of blank solution

Transfer an aliquot of dilution solvent in an amber vial and analyse by LC/MS.

5.5 Preparation of test solutions

Accurately weigh 75 mg (\pm 5 mg) of test sample in a 10 mL volumetric flask, dissolve with 5 mL of acetonitrile and dilute to volume with water.

Filter the slightly clouded solution through a 0.45 μ m PTFE or PVDF filter, transfer an aliquot of the clear solution in an amber vial and analyse by LC/MS.

The test solutions are stable for 3 days if stored at 20 °C and daylight in the autosampler tray.

6 Experimental
6.1 HPLC

| | |
|-------------------------------|---|
| HPLC column | Supelco Ascentis Express RP-Amide, 150 x 3.0 mm, 2.7 μ m Art.-no. 53919-U or equivalent |
| Column thermostat | 40 °C |
| Injection volume | 50 μ L |
| Injection mode | Injection with needle wash |
| Solvent needle wash | Water / acetonitrile (1:1) |
| Autosampler thermostat | 20 °C |
| Eluent A | 10 mM Ammonium acetate: Weigh 0.77 g \pm 0.01 g ammonium acetate in a 1 L bottle, add 1 L of water deionized and agitate until the ammonium acetate is dissolved. Afterwards adjust the pH with ammonia hydroxide solution to pH 8.0 – 8.5. |
| Eluent B | Acetonitrile |
| Flow rate | 0.5 mL/min |

| Time [Min] | A [%] | B [%] |
|-------------------|--------------|--------------|
| 0.0 | 55 | 45 |
| 3.0 | 55 | 45 |
| 8.0 | 5 | 95 |
| 13.0 | 5 | 95 |
| 13.1 | 55 | 45 |
| 16.0 | 55 | 45 |

6.2 Single quad mass spectrometer

| | |
|-------------------------|--|
| Ionization mode | API-ES |
| Polarity | negative |
| SIM Ion | 227.1 m/z Signal from 3.5 to 5.0 min From 0.0 to 3.5 min and 5.0 to 16.0 min → waste |
| Fragmentor | 170 V |
| Capillary voltage | 4500 V |
| Drying gas | 12 L/min |
| Drying gas temperature | 250 °C |
| Nebulizer pressure | 35 psig |
| Gain | 3.0 |
| Dwell time | 1000 msec |
| Expected retention time | 4.0 min |

6.3 Order of injections

The following order of injections is recommended:


| Sample description | Number of injections | Remarks |
|--|---|---------------------------|
| Blank solution, comparison solution respectively | n | System conditioning |
| Comparison solution 1 (CS1) LOQ | 1 | Linearity, LOQ |
| Comparison solution 2 (CS2) | 1 | Linearity |
| Comparison solution 3 (CS3) | 1 | Linearity |
| Comparison solution 4 (CS4) | 1 | Linearity |
| Comparison solution 5 (CS5) | 1 | Linearity |
| Comparison solution 6 (CS6) | 1 | Linearity |
| Blank solution | 1 | - |
| Test sample | n | - |
| Comparison solution 4 (CS4) | at least after every 10 th injection | SST – Performance control |
| Blank solution | 1 | - |

7 System Suitability Test

7.1 System Suitability Tests

| Test | Solution | Evaluation | Acceptance criteria |
|---------------------------|----------|---|---------------------|
| SST - LOQ | CS1 | S/N-ratio of Bisphenol A (SIM peak) | ≥ 10:1 |
| SST - Linearity | CS1-CS6 | Correlation coefficient (r ²) | ≥ 0.98 |
| SST – Performance control | CS4 | Recovery of Bisphenol A | 70 – 130% |

7.1.1 SST - LOQ - Signal to noise ratio



$$S/N = \frac{2 \cdot H}{h}$$

7.1.2 SST - Performance control

$$R_{\text{analyte}} = \frac{M_{\text{analyte}} \times 100\%}{T_{\text{analyte}}}$$

R_{analyte} Recovery analyte [%]

M_{analyte} Measured analyte content in comparison solution 1

T_{analyte} Theoretical analyte content in comparison solution 1

8 Evaluation and reporting


The concentration of Bisphenol A in the test sample is calculated with the external standard method (calibration function: linear).

This calculation is usually performed by the CDS.


The results will be reported in ppm (mg/kg).

9 Chromatograms


9.1 Blank solution



9.2 Comparison solution 1 (CS1) / LOQ / 0.5 ppm



9.3 Test solution (epoxy resin)



10 Abbreviations

| | |
|--------|--|
| API-ES | Atmospheric pressure ionization – Electrospray |
| HPLC | High performance liquid chromatography |
| MS | Mass spectrometer |
| m/z | Mass to charge ratio |
| CDS | Chromatography data system |
| CS | Comparison solution |
| RSD | Relative standard deviation |
| LOD | Limit of detection |
| LOQ | Limit of quantification |
| ppm | Parts per million (also mg/kg) |
| SIM | Selected (single) ion monitoring mode |
| S/N | Signal to noise |
| SOP | Standard operation procedure |
| SST | System suitability test |

11 Changes to previous edition

Changes

Validation completed and results considered

Chapter 5.3 and 5.5: Stability of comparison solutions and test solutions considered

Justification of changes

New edition after successful limit test validation

Impact on validation status

Limit test method is fully validated according to ICH Guideline Q2 (R1)